

ENDYMED

3DEEP SKIN SCIENCE



The unique and complete skin tightening, skin rejuvenation and body contouring system.

Technology

The technology behind ENDYMED® is second-to-none with the unique and patented 3DEEP Skin Science®, an advanced multi-source, phase-controlled radiofrequency technology.

3DEEP Skin Science allows the application of focused radiofrequency in different ways, treating the epidermis, the dermis or deep into the hypodermis with precision. The energy is controlled by the most advanced software, electronics and safety features available.

Clinical benefits

It's this industry-leading technology that enables the flexibility to treat every key skin rejuvenation indication, and deliver body contouring and tightening treatments - all with outstanding, clinically proven, long-lasting results and minimal pain or downtime for your clients:

- Lax and sagging skin on the face, including the peri-ocular area, cheek, jowls and neck
- Fine lines, deep lines and wrinkles
- Impaired skin texture, tone and luminosity

- Lax, sagging body skin, impaired body shape and cellulite on areas such as the tummy, flank, thighs and buttocks.

It is even possible to treat all of these indications from a single system - with a choice of platforms you can simply add treatments and handpieces over time.

ENDYMED TIGHTEN

Facial skin tightening and lifting

Powerful tightening and lifting of lax, sagging skin on the cheeks, jowls, neck and around the eyes^{1,2}. Using different handpieces to target specific areas, the radiofrequency is delivered deep into the skin to produce optimum results with no pain or downtime for your clients.

A course of 6 treatments is recommended and can be carried out by aestheticians NVQ3 and above.

ENDYMED V-SHAPE

Focused and deep skin tightening for the jawline and under the chin

The ultimate jowl busting treatment, ENDYMED V-SHAPE, delivers ultra-focused radiofrequency energy deep into the skin to target those stubborn areas of laxity – the jowls and submental area. It's the specialist MINI-SHAPER handpiece that focuses and intensifies the energy into the dermis without any pain or downtime for the client. It's ideal for use in combination with ENDYMED TIGHTEN or INTENSIF

A course of 6 treatments is recommended, and can be carried out by aestheticians NVQ3 and above.


ENDYMED FSR

Fractional skin resurfacing and tightening

ENDYMED FSR is the ideal treatment for fine lines and wrinkles, sun damaged skin and impaired texture and tone³. Here, 3DEEP Skin Science is used to rejuvenate the skin from the outside in, with controlled epidermal micro-resurfacing combined with volumetric dermal heating. Using the specialist FSR handpiece, treatments can be carried out for the full face or to target those common areas of concern – around the eyes and mouth. ENDYMED FSR can also be used to treat scars and stretch marks with fantastic results. The treatments are comfortable – with a topical anaesthetic applied beforehand – and there's minimal downtime for the client.

UNIQUE, ULTRA-FOCUSED
RADIOFREQUENCY TECHNOLOGY

ENDYMED 3DEEP
SKIN SCIENCE



A course of 3 treatments are recommended and carried out by advanced aestheticians or healthcare professionals.

ENDYMED INTENSIF

Ultra-deep fractional skin remodelling, tightening and lifting

ENDYMED INTENSIF rejuvenates the skin from the inside out, delivering the ultimate treatment for deep lines and skin laxity⁴, while remodelling the skin's underlying structure and restoring its vitality and radiance. During this clinically advanced treatment, radiofrequency is delivered three layers deep into the skin using unique, non-insulated gold needles and fractional pulse mode technology, creating fractional columns of coagulation and volumetric heating⁵. The specialist INTENSIF handpiece delivers comfortable treatments (with topical anaesthetic applied

beforehand), which remarkably have no downtime for the client.

A course of 3 treatments are recommended and can be carried out at healthcare professional led clinics.

ENDYMED SHAPER

Focused skin tightening and firming for the body

Ideal where skin laxity on the tummy, bottom, flanks or thighs are a concern, ENDYMED SHAPER delivers ultra-focused radiofrequency energy deep into the skin to firm, tighten and lift⁶ it. It's the SHAPER handpiece that focuses and intensifies the energy into the dermis without any pain or downtime for the client.

A course of 6-8 treatments is recommended and can be carried out by aestheticians NVQ3 and above.

ENDYMED CONTOUR

Ultra-deep body contouring and tightening with inch-loss, cellulite and fat reduction

The ultimate body contouring treatment, ENDYMED CONTOUR, combines high-power 3DEEP Skin Science with a sophisticated modulated vacuum massage. This powerful combination delivers radiofrequency energy ultra-deep into the skin, while the vacuum massage works synergistically to enhance the treatment of cellulite and fat reduction. The specialist CONTOUR handpiece is used to treat the tummy, bottom, flanks and thighs with no pain or downtime for the client.

A course of 6 treatments is recommended and can be carried out by aestheticians NVQ3 and above.

HOME USE-PROGRAMMES

For those wishing to enhance their clients' results and maximise ongoing retail sales, there are carefully designed, clinically proven skincare programmes for each of the ENDYMED facial rejuvenation treatments.

SEE OUR FANTASTIC CLINICAL RESULTS

ENDYMED TIGHTEN

Before



After



Photographs courtesy of I. Vider, M.D., Medical OR, Hertzlia, Israel

ENDYMED CONTOUR

Before



After



Photographs courtesy of Lilach Kaplan, Tel Aviv, Israel

REFERENCES

1. Harth, (2011). Journal of Cosmetic Dermatology.
2. Tanaka, (2013). Journal of Cosmetics, Dermatological Sciences and Applications, 3, 110-116.
3. Dahan et al, (2013). J Cosmet Laser Ther.
4. Tanaka, (2015). Lasers in Surgery and Medicine.
5. Harth, (2013) Data on File.
6. Royo Dela Torre, (2011). J Clin Aesthet Dermatol.

ENDYMED[®] and 3DEEP[®] are trademarks of ENDYMED[™] Medical Ltd, Israel.